

# Metal Tubes



## **Precellys Metal Tubes, for the homogenization of very hard samples and cryogrinding**

- Different tubes sizes
- Compatible with dry and liquid homogenization
- Suitable for Cryogrinding
- To be combined with Precellys beads

# METAL TUBES

The Precellys® Metal tube range ensures the homogenization of the hardest samples, while benefiting from the well-known efficiency provided by the 3D movement that defines the Precellys® instruments. Combined with the Precellys beads, they offer the highest performance.

The Precellys® Metal Tubes are also compatible with Cryogrinding, a technique that allows for the submerging and snap freezing of the tubes in liquid nitrogen.

## Cryolis® Evolution P000671-CLYS2-A.0



Discover our video on cryogrinding

## Beads



2.8mm Stainless steel beads (500g)	P000925-LYSK0-A.0
2.8mm Zirconium oxide beads (325g)	P000926-LYSK0-A.0
1.4mm Zirconium oxide beads (325g)	P000927-LYSK0-A.0
5,0mm Zirconium oxide beads (100 beads)	P000930-LYSK0-A.0
6,8mm Zirconium oxide beads (100 beads)	P000931-LYSK0-A.0

## 7mL Metal Tube + Holder P000951-LYSK0-A.0




Compatible with Precellys® Evolution and Precellys® Evolution Touch

## 2mL Metal Tube P002290-LYSK0-A.0



Compatible with Precellys® 24 Touch, Precellys® Evolution and Precellys® Evolution Touch up to 10 tubes simultaneously

## Discover our comprehensive range of solutions



Bertin Technologies SAS  
Parc d'activités du Pas du Lac  
10 bis, avenue Ampère  
78180 Montigny-le-Bretonneux, FRANCE  
01 39 30 60 00

WWW.BERTIN-TECHNOLOGIES.COM